



# Onsite Case Cart Repair

Herzog Surgical's onsite capabilities include case carts, O.R. tables, and I.V. poles and much more. Herzog Surgical's **SteelSure™** services will add life to your stainless steel equipment. **SteelSure™** will help save you from the hassle of replacing equipment as well as helping to meet budgetary goals.

## Repairable Items Include

- Autoclave Carts
- Autoclave Racks
- Case Carts
- Crash Carts
- IV Poles
- Surgery Tables
- Stainless Steel Items
- Items with Casters



## Repairs:

- Caster Replacement
- Welding
- Custom Fabricating
- Bumper Replacement
- Preventive Maintenance
- Doors
- Latch Assemblies
- Hinges
- Door Hinges, Racks & Shelf Brackets
- Cart Body (including dent removal, sticker removal)